



THE FUTURE IS WIRELESS

📍 135 W Central Rd. Suite 204
Schaumburg, IL 60195

☎ 224-801-5900

✉ Sales@A2ZImaging.com



www.Air2Zed.com

SMARTER DENTISTRY BEGINS WITH
SMARTER TECHNOLOGY

The Air2Zed system makes your existing X-ray sensors wireless. This plug-and-play solution expands the possibilities of your imaging workflow while requiring no modifications to your current sensors or software.

WIRELESS

Makes your imaging workflow wireless, efficient, and clutter free.



PORTABLE

Use a single sensor across multiple rooms.



EXPANDABLE

Designed to scale with your growing practice.



A WIRELESS IMAGING SOLUTION IS HERE!

- ✔ Plug-and-play
- ✔ 30-ft range
- ✔ No modifications to sensors
- ✔ Regain 40+ hours production time

Would Air2Zed benefit your office and improve overall efficiency?

- ✔ HYGIENISTS: 94%
- ✔ DENTISTS: 90%

*Sources: CDHA, independent unpaid surveys, national dental meetings

- ✔ No software to install
- ✔ Easily use one sensor across multiple rooms
- ✔ Supports most sensor brands
- ✔ Extend sensor life



THE FUTURE OF THE DENTAL OPERATORY.

The leader in dental technology innovation, Zinnov introduces Air2Zed, a sleekly designed wireless X-ray solution to keep your office safer, cleaner, and more efficient.

Say goodbye to tangled wires, expensive damage, marks on your floor, lack of mobility, patient discomfort, and overall office inefficiency

AIR2ZED EXPRESS STARTER SYSTEM

INCLUDES A SYSTEM TO ACCOMMODATE 1 SENSOR AND 2 COMPUTERS/OPS

Supports one sensor and up to two workstations



• One transmitter



• Two receivers



• Charger and accessories included

HOW IT WORKS

MAKE YOUR EXISTING X-RAYS WIRELESS WITH AIR2ZED IN **3 SIMPLE STEPS.**

There is no additional software to install and no modifications required to your sensors. It truly is a plug and play system.

STEP 1:

Connect a Zed(receiver) to a workstation.



STEP 2:

Connect your X-ray sensor to the portable transmitter.



STEP 3:

Click the corresponding button to pair with the room of your choice



THAT'S IT!

Your x-ray sensor is now wirelessly paired with the workstation, and you can continue to follow your normal workflow.

And the best part is, the Air2Zed system is designed to grow with your practice by adding pieces as needed!

HOW TO BUILD YOUR AIR2ZED

01



Begin by selecting an Air2Zed Express Starter System – This will support 1 sensor and up to 2 computers.

02

Choose how many ADDITIONAL sensors you would like to add to your system. Add 1 Air2Zed Express Transmitter for each of these additional sensors.

(Example; If you currently have 2 sensors that you use, you would need 1 additional transmitter, as the first is included in the starter system.)



03



Choose how many ADDITIONAL computers/monitors you currently send images to (Generally, offices use one in each operator). Add 1 Air2Zed Receiver (Zed) each of these additional computers.

(Example: If you are sending images to 5 computers, you would need 3 additional receivers, as the first 2 are included in the starter system)

In the above scenario, if you use 2 sensors in 5 ops, you would purchase 1 Air2Zed Express Starter System, 1 Air2Zed Express Transmitter, and 3 Receivers (Zeds).

SUPPORTED SENSORS

Visit website for most current list



vatech



Carestream



PARTNERS

Visit website for most current list



FAQs

HOW MANY IMAGES CAN BE TAKEN BEFORE THE BATTERY RUNS OUT?

The battery in your Air2Zed is designed to easily last through your entire day. When you close for the day, simply set it to charge and it will be ready to go again when you open in the morning.

CAN PIECES BE BOUGHT INDIVIDUALLY?

Yes, they can. Most customers begin with our starter system and then can expand the system to their needs, as little as one piece at a time.

DO THE ZEDS (RECEIVERS) NEED A BATTERY?

Because your Zeds are powered by the USB port, they do not require a battery.

WHAT IS THE WARRANTY?

Air2Zed comes with a standard 2 year warranty. Visit website for details.

HOW LONG DOES IT TAKE TO COMPLETELY CHARGE THE BATTERY?

The Air2Zed transmitter has a quick charge battery that can be fully recharged from empty in approximately 2 hours.

HOW SHOULD THE DEVICE BE DISINFECTED?

Clean your Air2Zed with the same sanitary wipes that you regularly use in your office. Air2Zed should not be submerged in water or other liquid.

WILL THE DEVICE WORK WHILE WEARING LATEX GLOVES?

Yes. All buttons have a textured feel, making it simple to operate.

IS THE DATA MODIFIED?

No, the data from the X-ray sensor is not modified, compressed or altered in any way.

Enhance your x-ray procedure with our cutting-edge wireless technology, Air2Zed. Request a personalized product demonstration today to experience the game-changing benefits of making YOUR CURRENT x-ray sensors wireless, firsthand. Reach out to your trusted equipment dealer and ask for Air2Zed. Discover how Air2Zed can streamline your imaging workflow, improve efficiency, and elevate patient care. Don't wait to unlock the potential of your favorite sensor!

www.Air2Zed.com

POWERED BY ZINNOVI



*Proudly designed and
manufactured in the USA*



121 Exchange Blvd, Glendale
Heights, IL 60139



888-629-6646



www.dentalsalesinc.com

